

Covetrus Dental Brochure 2024

Browse through our selection of dental products available for both small animal and equine dentistry.



Contents

Small Animal

Small Animal

MOUTH GAGS & THROAT PACKS	3
SCALE AND POLISH	3
EXTRACTION	5
BURS	6
HANDPIECES & ACCESSORIES	6
Home Care	
DIET	7
TOOTHBRUSHING & PASTES	8

Equine



MOUTH GAGS & LIGHTS	11
DRENCHING SYRINGES	11
PICKS & PROBES	11
MIRRORS	12
FORCEPS	12
FLOAT HANDLES, BLADES & FILES	12



MOUTH GAGS & THROAT PACKS

- Easy to use and clean
- Plastic/rubber inserts
- Spring
- Self-locking

Limited stock



2801183
CVET Mouth Retractor For Cats 6cm


2801184
CVET Mouth Retractor For Dogs Large 15cm

2801186
CVET Mouth Retractor For Dogs Small 7cm

2801172
CVET Jaw Retractor 13cm

2801240
CVET Rabbit Jaw Retractor
Specifically designed for rabbits.

Limited stock



ROTHR04
Robinson Dental Throat Pack Small 5.1 x 5.1cm 80pk

ROTHR05
Robinson Dental Throat Pack Medium 6.3 x 6.3cm 50pk

ROTHR06
Instrapac Throat Pack Medium 2's 40pk

Absorbent cotton pad with securely attached removal cord for safety.



SCALE AND POLISH

2801831
CVET Vet Piezo Scaler LED
Compact ultrasonic scaler unit specifically designed for veterinary dentistry.

- Scaling, endo and automatic water supply functions
- With adjustable intensity and water volume
- Detachable autoclavable handpiece and tip change
- Easy to clean



KRUUSE Piezo Scaler is supplied with:

- 271400, 271401, 271402 3 Bonart scaling tips No 1, 2 & 3
- 390536 Bonart Tip wrench
- 271190 Pump, 5 l water bottle

- 271103 Foot pedal control, on and off
- 271395 Mikromotor E fitting and holder
- 282464 Straight nose cone 1:1

- 282485 KRUUSE Starter Set for polishing (KRUUSE Oscillating Prophy Angles, Polishing Paste, Polishing Cups)



KRSCA10
KRUUSE ART-SP-2 Piezo Scaler & Micromotor



KRSCA09
TIP BS-1 for KRUUSE Piezo Scaler



KRSCA16
Tip BS-3 for KRUUSE Piezo Scaler Satelec P3II



KRSCA15
Tip #1 Mk1 Piezo Scaler BE1 (SP2 & P3, EMS)


KRSCA11
NSK Scaler Tip G 4

KRSCA12
NSK Scaler Tip G 2


KRSCA13
NSK Scaler Tip G 3

KRSCA14
NSK Scaler Tip G 6

NSK scaler tip G 4	• Removal of supra gingival tartar. • Removal of marginal gingival calculus.			
NSK scaler tip G 2	• Removal of supra gingival tartar.			
NSK scaler tip G 3	• Removal of supra gingival tartar.			
NSK scaler tip G 6	• Removal of supra gingival tartar. • Removal of tartar narrow interdental space. • Removal of marginal gingival calculus • Removal of gingival biofilm under natural tooth			



CANINE DENTAL CHART



FELINE DENTAL CHART

2801250
CVET Scaler Jaquette 16.5cm
For the removal of supragingival calculus and plaque (above the gum line).



2801076
CVET Dental Explorer #17/23 Double-Ended Comfort
Double-ended explorer used to test tooth stability and the condition of the gum margins.



The CVET 'Comfort' instruments are designed to fit well in the hand to prevent discomfort during repetitive hand movements.

Curettes for the removal of deep subgingival calculus and for the finishing of subgingival root surfaces.

2801071
CVET Periodontal Curette Gracey #1-2 Comfort
2801073
CVET Periodontal Curette Gracey #3-4 Comfort



KRDEN59
KRUUSE Oscillating Prophy Angle With Cup
2800791
CVET Disposable Prophy Angles Universal 100pk

- Eliminates splatter of prophy paste
- Avoids cross-contamination
- Eliminates frictional heat, making it friendlier on the surface of the tooth while protecting enamel



2800819
CVET Prophy Paste Cups Medium Grit 200pk
Contains a unique blend of polishing and cleaning agents. Fluoride free. Non-splatter formula.

IMDEN09
Dental Scaler Paste Pink 637820 10pk
MAPRO03
Prophy Paste Fine Grit 5pk

2800167
CVET Tongue Depressors Non Sterile 15cm 500pk
Used to examine mouth and throat.



2800862
CVET Plaque Indicating Swabs 25pk
Useful when scaling teeth to ensure that all the plaque has been removed during the procedure.



AASAN01
Sanos Vet Dental Sealant Single
Gingival dental sealant for dogs and cats that seals the subgingival gum line against the formation of plaque. Clinically proven to last for 6 months without any take home application required.



2800812
CVET Pet Fluoride Foam Vanilla 125g
Fluoride solutions are antibacterial and can reduce dental hypersensitivity on exposed root surfaces.



EXTRACTION

2801235
CVET Periosteal Elevator Molt #9 171mm Comfort



Used to lift full thickness soft tissue flaps.

2801234
CVET Periosteal Dental Elevator #A2 171mm Comfort

2801085
CVET Elevator Bein 2mm

2801086
CVET Elevator Bein 3mm

2801087
CVET Elevator Bein 4mm




Used to loosen the tooth from the periodontal ligament and ease extraction.

KRRET03
KRUUSE Minnesota Retractor
Used to retract and protect soft tissue when working with drills.




2801247
CVET Root Elevator Hylin 3mm
Used for the extraction of a tooth or retained roots.



2801246
CVET Bone Rongeur Forceps Luer Slim Curved 15cm

2801065
CVET Bone Rongeur Forceps Mini Friedman 12.5cm



Used for the removal of bone following tooth extraction.

2801112
CVET Extraction Forceps Large Breed

2801113
CVET Extraction Forceps Small Breed



For clean and precise tooth extractions.

KREXT07
KRUUSE Extraction Forceps Deciduous Teeth



Luxators are fine sharp tools designed to remove the carnassial without the need to split the tooth.



VNLUX01
Luxator Svenska 1mm for Feline Teeth
HSLUX02
Sdi Luxator Straight 5mm
HSLUX01
Sdi Luxator Straight 3mm

2801177
CVET Luxator for Rodent Molars



2800079
CVET Small Surgery Set Sterile

- 2 x Surgical drape 45 x 75 cm
- 1 x Table cover 75 x 90 cm
- 1 x Dental drape 160 x 150 cm
- 2 x Cord drape for tubing 8 x 120 cm
- 2 x OP adhesive tape 10 x 50 cm
- 1 x Wrapping sheet 75 x 80 cm

BURS

MABUR56
Round Tip Bur #2 FG 19mm
Carbide 5pk

MABUR57
Round Tip #2 FG Surgical 25mm
Carbide 5pk

MABUR32
Round Tip #6 FG Surgical Bur 25mm
Carbide 5pk

MABUR33
Round Tip 1/4 FG Surgical Bur 25mm
Carbide 5pk

MABUR34
Round Tip #4 FG Bur 19mm
Diamond 5pk

MABUR39
Round Tip #8 FG Bur 19mm
Diamond 6pk

MABUR41
Round Tip #4 FG Bur 19mm
Carbide 5pk

MABUR42
Round Tip #4 FG Surgical 25mm
Carbide 5pk

MABUR43
Round Tip #1 FG 19mm
Carbide 5pk

MABUR60
Round Tip Bur #6 FG 19mm
Carbide 5pk

MABUR45
Tapered Fissure Bur #701 FG 19mm
Carbide 5pk

MABUR46
Tapered Fissure Bur #701L FG 19mm
Carbide 5pk

MABUR48
Tapered Fissure Bur #701L Surg.FG 25mm
Carbide 5pk

MABUR50
Tapered Fissure Bur #699 FG 19mm
Carbide 5pk

MABUR51
Tapered Fissure Bur #700 FG 19mm
Carbide 5pk

MABUR52
Tapered Fissure Bur 702 Surg.FG 25mm
Carbide 5pk

MABUR53
Tapered Fissure Bur #702 FG 19mm
Carbide 5pk

MABUR55
Pear-Shape Bur #330 FG 19mm
Carbide 5pk

MABUR58
Pear-Shape Bur #332L FG 19mm
Carbide 5pk

MABUR59
Pear-Shape Bur #331 FG 19mm
Carbide 5pk

WHICH BUR TO USE?

ROUND: used for cavity preparation, creating access points undercuts and channels for luxator blades in extraction.

TAPERED FISSURE: used for sectioning multi-rooted teeth and reducing crown height.

PEAR: Used for cavity preparation, access points and splitting roots of small teeth.

KRBR07
KRUUSE FG Bur Set
445mm 10pk



HANDPIECES & ACCESSORIES

Top quality handpieces from BA International and KRUUSE.



2801568
BA Contra angle
1:1 Handpiece-latch
type bur release



2801571
BA BA45E Contra
angle 1:1 Handpiece
push button/bur



2801576
BA BA175E Straight
1:1 Handpiece.
No Spray



2801570
BA Prophy Contra-
angle 1:1 for snap
on cups/brushes



KRDEN05
Metal Prophy angle
271127



KRDEN55
KRUUSE
High-Speed
Handpiece
W/LED

KRUUSE High-Speed handpiece has automatically generated LED light run by the air in the unit.



KRDEN56
NSK Low-Speed Motor
282462



KRDEN58
NSK Straight Nose Cone
1:1 282464

Scan for tips
on instrument
cleaning and
maintenance



KRDEN14
KRUUSE Dental Formula
250ml 6pk
Chlorhexidine dental solutions
serve as a coolant for the
moving parts of turbines.
Also reduce bur cutting time
and help maintain cleanliness
in the mouth.



2801497
CVET Dental Solution
250ml

2801492
CVET Oil (Clipper) Spray
With Nozzle 500ml
Universal oil spray for straight
and contra-angle handpieces.
Degreases, cleans, and
lubricates. Includes spray
head.



KROILO1
NSK Oil Bottle 5ml
Small oil drip bottle for precise
application. It is advisable to
apply oil every time before
autoclaving any handpiece
or air-driven motor.



MALUB07
Bio Lube Lubricant
Pen 25oz
Synthetic
biodegradable
handpiece lubricant
that is safe for all high
speed and low speed
handpieces.



2801306
CVET Sharpening
Stone Arkansas
Fine Grit 15x5cm
Restores the
cutting edge on
dull instruments.
Fine grit,
rectangular.



Home Care

DIET

2802249
CVET NutriCareVet Dental
Support Soft Chew 60pk
CVET dental soft chews
provide support for
healthy teeth, gums and fresh
breath. Developed for dogs
and puppies from 12 weeks.



2800745
CVET Petcare
Ostrich Twister
Chew 10pk
100% all-natural dog
chews made from
the tendons of South
African Ostriches.



NVTRE74
Calibra Joy Dog Classic
Dental Brushes 85g
A healthy and easily
digestible treat helps
fight dental plaque and
calculus.



Extremely palatable
rawhide and collagen
based chews.

CELOG25
Logic Prozym S 15 Chew Strip
CELOG26
Logic Prozym M 15 Chew Strips

Delicious supplements for dogs with added benefits, catering to various tastes and sizes. They effectively control plaque and tartar, promoting fresher breath, while also providing omega-3 and -6 fatty acids and antioxidants for immune support.



MRPLA20
PlaqueOff
Dental Cat
Bites 60g



MRPLA12
PlaqueOff Animal
Dental Bites
60g

MRPLA13
PlaqueOff Animal
Dental Bites
150g



MRPLA16
PlaqueOff Dental Bones
Veg Fusion 482g

MRPLA17
PlaqueOff Dental Bones
Turkey & Cranberry 485g

MRPLA18
PlaqueOff Dental Bones
Natural Bacon 482g

MRPLA19
PlaqueOff Dental Bones
Chicken & Pumpkin 485g

MRPLA21
PlaqueOff Mini Bones
Peanut Butter & Banana 340g



TOOTHBRUSHING & PASTES



Prevention is key

Maintaining excellent oral health relies on daily oral hygiene practices, with brushing considered gold standard. Additionally, daily chewing activities can contribute effectively to oral health maintenance.

Take a look at what we recommend for pet parents:

Finger toothbrushes are the easiest type of brush for beginners.



2800664 CVET Dental Care Finger Toothbrush 2pk

With zinc gluconate.



2800817 CVET Pet Toothpaste for Animals 60g

Hexarinse is an oral rinsing solution which helps to control levels of bacteria in the mouth.



SJHEX01 Virbac Hexarinse Chlorhexidine 237ml

Drinking water additive. Innovative formul which inhibits plaque forming bacteria.



VIAQU04 Vet Aquadent 250ml
VIAQU03 Vet Aquadent 500ml

Helps remove and prevent dental plaque in cats and dogs. Use with or without brushing.

ACORO01 Orozyme Oral Hygiene Gel 70g
ACORO02 Orozyme Oral Hygiene Gel 70g 12pk



2800815 CVET Pet Chlorhexidine 0.12% Cleaning Gel 114g
Chlorhexidine products help to reduce the number of oral bacteria that causes both periodontal disease and bad breath in cats and dogs.



VJSTO01 Stomodine F Gel

Antibacterial / soothing gel to help maintain healthy teeth and gums in dogs and cats. Contains Chlorhexidine and Tris EDTA to help manage bacteria that cause plaque build-up and bad breath. Phytosphingosine and Safflower Oil combine to soothe the gums and quickly reduce redness.



Based on a multi-enzyme patented system, Logic Oral Hygiene Gel supplements the animal's own defence mechanisms to help fight harmful bacteria in the mouth.

CELOG04 Logic Oral Hygiene Gel 70g Single
CELOG11 Logic Oral Hygiene Gel 12pk Display Pack



SJTOO31 Enzymatic Toothpaste Poultry Dog 70g

Enzymatic Toothpaste for Cats Fish Flavour 43g
Contain CET Dual Enzyme, an innovative formula which inhibits plaque forming bacteria.

SJTOO34 Enzymatic Toothpaste Kit Dog Poultry

SJTOO35 Enzymatic Toothpaste Kit Cat Fish
Include 1 x toothpaste, 1 x finger brush and 1 x dual-ended brush



FODEN01 Dentisept Paste Dog & Cat 20g Single
Vanilla-flavoured palatable oral paste containing chlorhexidine.

Plaqtiv+® is a new oral care range for cats and dogs that delivers real-world results with its patented Xpersiv™ technology. Plaqtiv+® captures minerals, inhibiting biofilm maturation and thus formulation of calculus.



ACPLA02 Plaqtiv+ Oral Care Toothpaste Malt 70g



ACPLA01 Plaqtiv+ Oral Care Water Additive Original 500ml



ACPLA07 Plaqtiv+ Oral Care Oral Spray Vanilla Mint 60ml



ACPLA08 Plaqtiv+ Oral Care Hygiene Wipes Vanilla

ProDen PlaqueOff® Powder is your 100% natural and easy solution to pet oral problems – just sprinkle the powder once a day on wet or dry food to help reduce problems and pain in the future, as well as freshen bad breath.



MRPLA01 PlaqueOff Animal 60g
MRPLA04 PlaqueOff Animal 180g
MRPLA02 PlaqueOff Animal 420g



MRPLA11 PlaqueOff Cat 40g



DNTOO01 Roast Dinner Toothpaste



Kits contain:

- Petosan Double-Headed Pet Toothbrush
- Petosan Toothpaste with Poultry Flavor
- Microfibre Finger brush

- PST0012**
Petosan Dental Kit Toy Dog/Puppy
- PST0009**
Petosan Dental Kit Small
- PST0008**
Petosan Dental Kit Medium
- PST0007**
Petosan Dental Kit Large



PST0002
Petosan Microfibre Tooth Cleaner

This antimicrobial, chemical-free, self-cleaning microfibre finger cloth easily slips over your index finger and has a strap to attach to your middle finger ensuring easy cleaning of teeth and mouth.



The brush heads are angled so that the bristles meets the tooth at a 45-degree angle. This is ideal in order to be able to clean at and below the gum line where the harmful bacteria thrive.

- PST0011**
Petosan Toothbrush Toy Dog/Puppy
- PST0005**
Petosan Toothbrush Small Dogs
- PST0004**
Petosan Toothbrush Medium Dogs
- PST0003**
Petosan Toothbrush Large Dogs



PST0001
Petosan Poultry Flavour Toothpaste 70g
The Petosan Pet Toothpaste is designed for use by those who brush their dog's teeth every day. It is formulated to be very tasty, helping to make brushing a positive experience.



The Arm & Hammer oral care for dogs range provides an effective and convenient way to care for your pet's teeth and gums. Each product is specifically designed to cater for the requirements of individual dogs and their owners.



- DMDEN11**
Arm & Hammer Coconut Dental Kit For Dogs
- DMDEN08**
Arm & Hammer 360 Degree Toothbrush Large
- DMDEN09**
Arm & Hammer 360 Degree Toothbrush Small
- DMDEN14**
Arm & Hammer Coconut Mint Toothpaste For Dogs
- DMDEN13**
Arm & Hammer Coconut Mint Toothpaste For Puppies
- DMDEN10**
Arm & Hammer Coconut Water Dental Spray 4oz
- DMDEN15**
Arm & Hammer Coconut Water Additive 16oz Dogs
- DMDEN16**
Arm & Hammer Coconut Water Additive 16oz Puppies

Equine

Scan the QR code for more on Equine Dentistry from KRUISE



MOUTH GAGS & LIGHTS

KREQU65
EQUIVET Millenium Speculum
The top quality stainless steel mouth speculum gives adequate access to the oral cavity. Easy to adjust through 15 settings with front actuated release levers. Comes with 2 kinds of biteplates (standard and leather). Biothane straps are water proof and flexible. They are easy to use and clean.



KRSPE06
EQUIVET Climax Speculum Gag
This high quality speculum is manufactured in forged stainless steel with adjustable biothane straps. Spring loaded ratchets with five steps to keep the mouth opened at the desired setting during use. Includes one set of standard mouth plates and a set of leather sleeved bar plates.



KRGAG07
Comet Mouth Gag / Wedge Plastic



KRRET03
KRUISE Minnesota Retractor
Used to retract and protect the tissue and the user when working with a drill. Each end also has a specific function: One holds back the cheek, while the other holds down the tongue.



KREQU82
EQUIVET Extra Wide Bit Plate (Pair)
These bit plates create additional working space by widening the side parts, providing a significant advantage when dealing with geriatric or overjet horses.



KRDEN06
EQUIVET DentLite
Developed for veterinary surgeons and equine dental technicians who need an efficient working light when operating in the oral cavity of the horse.

- High light intensity
- Correct colour representation
- No annoying cable between light source and battery



DRENCHING SYRINGES

For dosing large quantities of water to clean the horse's mouth before oral examination.

KRSYR37
KRUISE Mouth Flushing Dosing Syringe 300ml



KRSYR02
Mouth Wash Syringe Long Nozzle 400ml



PICKS & PROBES

KRDEN13
EQUIVET Dental Probe Fine Tip
This sturdy, long handled, fine tipped calibrated probe is ideal for evaluating the pulp horns and infundibula on the occlusal surface of cheek teeth.


KRDEN08
EQUIVET Dental Probe Heavy Tip
This long, heavy tipped calibrated probe is engineered to probe, clean, and measure the depth of periodontal pockets involving the equine cheek teeth.

KREQU43
EQUIVET Dental Irrigation Pick
This specially designed dental pick incorporates water irrigation to expedite the process of cleaning interdental spaces and periodontal pockets.

Considered the best choice for treating diastema.



KRDEN07
EQUIVET Dental Pick Set
• All heads can be rotated.
• The chisel head can be used to separate the gingiva from the sides of the tooth prior to extraction.
• The long and short osteotome heads can be used to separate teeth or tooth fragments.
• The long and short elevator heads can be used in the short hand piece to aid in removal of wolf teeth or they can be used on the long handle to aid in the removal of dental root fragments from the alveolus.
• The loop curette is an ideal tool to snag loose tooth or root fragments from the alveolus.
• Chuck with short and long pins is used as dental pick for probing alveolar pockets and areas of infundibular decay as well as for removing feed and debris.



MIRRORS



KRMIR10
EQUIVET Dental Mirror Large Model 50cm 240744
 A necessity to perform proper oral cavity examination. The mirror is a 50 cm long T action handle. The mirror is 48 mm. The angle is 30°.

KREQU63
EQUIVET Spare Mirror For The EQUIVET Dental Mirror
 Spare mirror for EQUIVET Dental Mirror with adhesive.



2801351
CVET Tooth-Mirror Fahrenkrug Horses & Large Animal
 Essential for oral examination. Tooth-mirror with handle for larger animals. Stainless steel.

FORCEPS

KREQU77
EQUIVET Diastema Periodontal Forceps

- Sturdy
- Long-handled
- Right angled
- 41 cm



KRFOR17
Reynolds Cap Extractor Forceps Lower Jaw 37.5cm
 Instrument for retained deciduous teeth (caps) in horses. Also works for small premolars and incisors.



FLOAT HANDLES, BLADES & FILES

KRBUC02
EQUIVET Dental Rasp Bucket Buddy - holds 8
 This is designed to protect the rasps while they are placed in the bucket. The protector consists of 8 tubes and it still gives access to the cleaning brush in the bucket.



KRBLA36
EQUIVET Carbide Blade Fine cut 19 x 55mm
 The EQUIVET MagFloat Carbide Blades have a design where the carbide part is very slim and glued to a highly magnetic plate.



KRBLA37
EQUIVET Carbide Blade Medium cut 19 x 55mm
 The carbide part is very slim and glued to a highly magnetic plate. Cannot be re-sharpened.



I have been using the MagFloat set in my practice and for undergraduate teaching since receiving them from KRUUSE and DLC Australia. They are an exceptional set of hand dental instruments which are designed to work effectively in all areas of the mouth. It is such a simple procedure to change the float blades and the lightweight instruments reduce the muscular fatigue on the practitioner even after treatment of multiple horses. These instruments have become my instruments of choice whenever hand instrumentation is required.

Gary Wilson BVSc MVSc MACVSc CMAVA
 Cert Teach, Associate Professor Queensland, AUS



EQUIVET MagFloats are designed so that changing blades is a very easy procedure as the blades are held on to the head of the float with very strong magnets. The blade is easily pushed away from the magnets with the T-blade punch when changing to a new blade or switching from a "push" to a "pull" function.

- Lightweight
- Reduces physical strain
- Rasp plate is held in place by very strong magnets
- Easy to change blades

KRRAS31
EQUIVET Mag Float Upper Premolar (6-8)
 Length end to end - 35.5cm

KRRAS33
EQUIVET Mag Float Lower Cheek (6-11)
 Length end to end 53cm

KRRAS35
EQUIVET Mag Float Upper & Lower Molar (11)
 Length end to end 61cm



KRRAS37
EQUIVET Float Dental Rasp 47cm
 Stainless steel with tungsten carbide chips at both ends and only on one side to avoid injury to the buccal mucosa and tongue.



KRBUC01
EQUIVET Stainless Steel Bucket 10 Litre with Brush
 10 litre stainless steel bucket. Supreme quality, specially designed for use with equine dental treatment.

KRBRU05
EQUIVET Nylon Brush For Dental Bucket
 For use with stainless steel bucket.